

Berlin, 17.04.2025

## **Cleaning and Disinfection of Colposcopes**

Dear Customer,

Please pay attention to the recommendations below.

### **Cleaning**

Routinely wipe the colposcope and the stand with a damp (not wet) cloth. Solutions in the usual dilution can be used for cleaning. Then wipe with a dry, lint-free cloth.

Cleaning solutions can be:

- Universal cleaner
- Hygiene cleaner
- Surface cleaner

### **Disinfection**

Disinfect the unit with conventional disinfectants in accordance with national hygiene standards and in accordance with your internal hygiene plan.

Conventional disinfecting agents can be:

- Ready-to-use quick disinfection preparation based on alcohols, eg. 2-propanol (isopropanol);
- Disinfectants based on amines, such as chloramine-T;
- Disinfectants based on quats, eg dimethyl didecyl;
- Active oxygen-based disinfectant (percompounds), eg. oxygen releasers.

The cleaning solutions and disinfectants should be used according to the manufacturers' instructions in consideration of the recommended exposure times.

---

**LEISEGANG Feinmechanik-Optik GmbH**

To avoid damage to the equipment, observe the following instructions for cleaning:

- To clean the lenses you can use an aqueous soap solution. To dry the lenses you should use cloths for optical equipment.
- Please never use abrasive or aggressive cleaning and disinfectant agents as this may cause damage and corrosion to the surface.

The power supply is not protected against the penetration of liquids. Penetrating liquid may cause electrical shock.

- Avoid drops and splash water as this may cause damage to the device.
- Do not spray the power supply with cleaning or disinfecting agents.



**Falko Sandmann**  
Quality Manager  
LEISEGANG  
Feinmechanik - Optik  
GmbH  
Leibnizstraße 32  
D - 10625 Berlin  
Tel.: (030) 31 90 09 - (0)  
Fax: (030) 31 90 92  
Leisegang Feinmechanik-Optik GmbH  
A Cooper Surgical Company